

Quality & Calibration Report



131 Business Center Drive, A-3
Ormond Beach, FL 32174
TEL: 877.678.4259 | 386.310.4933 M-F 9-6pm EST
FAX: 866.422.2356
EMAIL: sales@co2meter.com | www.co2meter.com | www.metroco2.com

Certificate of Calibration and Testing

Date: [REDACTED]
Order Number: 5017
Customer: [REDACTED]

Item	Qty	Serial Number	Calibration Gas	Reading	Data-Logging Test
CM-0003 SAMPLING DATALOGGER 100% CO2	1	109026	Nitrogen 100% (± 0.015%) P/N 221L-114 LOT # HAP-114-14	0 ppm (± 10 ppm)	Pass
			CO2 95% (± 5%) P/N 221L-35-95 LOT # CAM-35-95-1	95% CO2 (± 5%)	Pass

QA Inspector

We are committed to providing a quality product.

This instrument has undergone rigorous testing throughout its manufacture.

This instrument was assembled in USA, and has been factory inspected, Tested & Calibrated in accordance with the conditions and requirements of our registered Quality System, Operating Standards & Sales Agreements.

All calibration gases are traceable to industries standards.

NOTE : This certification is valid for 180 days from calibration date.

Calibration of this device is traceable to NIST